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Tel 0039-0429-788256 Fax 0039-0429-788560

E mail fabio.manoni@ulss17.it

CHEMICAL, MORPHOLOGICAL AND MICROBIOLOGICAL URINE EXAMINATION: PROPOSAL OF GUIDELINES FOR PRE ANALYTICAL STANDARDIZATION.

Summary:

These guidelines for the definition of preanalytical phase of urinalysis have been drafted under the auspices of the Italian Society of Laboratory Medicine (SIMeL) and of the Italian Society of Clinical Biochemistry and Clinic Molecular Biology (SIBioC).

Urinalysis, including both morphological assay and microbiological examination should always be performed on the basis of an appropriate medical request and the relative analyses should be selected according to the clinical needs, the patients selection, and the available technology as well. The reference materials for comprehensive urine examinations should include clear and intelligible recommendations for the preexamination, examination and postexamination processes. A correct preexamination process should include accurate information about the preparation of the patient, collection, handling, preservation, storage and transport of the specimen. A description of the local skinless and therefore the hierarchical level of urinalysis is also recommended.

The chemical-morphological urine exam has undergone some changes over the last few years. As such, it might be the time to update and modify accordingly the various activities throughout the total testing process. The advent of new technologies, which have made the laboratory report much more significant as regards the morphological features of this test, have also promoted some reflections about the necessity of further refining and improving the quality of the preexamination process, on the opportunity of consolidating and standardizing the examination process, but also on redefining the clinical objectives by developing a complete, integrated and more clinically meaningful report.

These guidelines have been elaborated and drafted on the basis of the current literature as well as on the personal experience of the authors.

The more significant aspects of these guidelines concerne the following aspects:

- Urinalysis should be performed on the basis of a specific clinical need
- The sample recommended for urinalysis is the first morning urine collected by a voided midstream technique
- Each Service of Laboratory Medicine should make available to the stakeholders printed instructions (with illustrations, whenever possible) concerning the appropriate procedure for collection of the specimen
- The use of standard devices for urine collection is recommended
- The use of an appropriate labelling system is recommended
- Urinalysis should be performed as soon as possible after collection. If examination is to be delayed by more than four hours, then the samples should be refrigerated
- The preparation of adequate protocols for handling, transport and storage of the samples is mandatory
- It is necessary to transmit preexamination information from the collection sites to the Laboratory that performs the analysis. The occurrence of violations from the standard should always be recorded.
- The hierarchical analytical level available in each laboratory should be declared

We hope that these guidelines concerning the preexamination process of urinalysis have a positive impact on the validation of local working practice, producing a significant improvement of the total quality of this test (and in the whole field of laboratory medicine as well), so that the results of testing will reflect the real condition of the urinary system of the patient.

Key words: Guideline, urinalysis, preexamination process.

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- ECLM- *European Urinalysis Guidelines*. The Scandinavian Journal of Clinical and Laboratory Investigation. 60: 1-96, suppl. 231, 2000⁴.
- CLSI GP-16 A3 Urinalysis and Collection, transportation, and Preservation of Urine Specimens; Approved Guideline – third Edition vol.29: n4: 4-21⁵.
- Scottish Intercollegiate Guidelines Network (SIGN). Management of suspected bacterial urinary tract infection in adults. A National Clinical Guideline. Edinburgh (Scotland), Scottish Intercollegiate Guidelines Network (SIGN); July 2006. SIGN publication n. 88⁶.

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Bibliografia

1. Coomarasamy A. Searching for evidence to inform clinical practice. *Current Obstetrics & Gynaecology* 2004; 14:142-6.
2. Lilford R, Richardson A, Stevens A, Fitzpatrick R, Edwards S, Rock F, et al. Issues in methodological research: perspectives from researchers and commissioners. *Health Technol Assess* 2001; 15:1-57.
3. Grilli R. AGREE uno strumento per la valutazione della qualità delle linee guida. Dossier 60 Bologna. Agenzia Sanitaria Regionale dell'Emilia-Romagna 2002.
4. ECLM - *European Urinalysis Guidelines*. *Scand J Clin Lab Invest* 2000; 60:1-96.
5. CLSI GP-16 A3 Urinalysis and Collection, transportation, and Preservation of Urine Specimens; Approved Guideline – third Edition vol.29; n4:4-21.
6. Scottish Intercollegiate Guidelines Network (SIGN). Management of suspected bacterial urinary tract infection in adults. A National Clinical Guideline. Edinburgh (Scotland), Scottish Intercollegiate Guidelines Network (SIGN); July 2006. SIGN publication n. 88.
7. Atkins D, Best D, Briss P, Eccles M, Falck-Ytter Y, Flottorp S, et al. Grading quality of evidence and strength of recommendations. *BMJ* 2004; 328:1490.
8. Guyatt G, Gutterman D, Baumann MH, Addrizzo-Harris D, Hylek EM, Phillips B, et al. Grading strength of recommendations and quality of evidence in clinical guidelines: report from American College of chest physicians task force. *Chest* 2006; 129:174-81.
9. Kouri T, Gyory A, Rowan M. ISLH recommended procedure for the enumeration of particles in urine. *Lab Hemat* 2003; 9:58-63.
10. Aspewall O, Hallander H, Gant V, Kouri T. European guidelines: a collaborative document produced by European clinical microbiologists and clinical chemists under ECLM in collaboration with ESCMID. *Clin Microb Infect* 2001; 7:173-8.
11. British Columbia Guidelines and Advisory Committee. Protocol for macroscopic and microscopic urinalysis and investigation of urinary tract infections. 2004. www.healthservices.gov.bc.ca/msp/protoguides.
12. Fogazzi G, Garigali G, Pirovano B, Muratore M, Raimondi S, Berti S. How to improve the teaching in urine microscopy. *Clin Chem Lab Med* 2007; 45:407-12.
13. Kouri T, Malminiemi O, Penders J, Pelkonen V, Vuotari L, Delanghe J. Limits of preservation of samples for urine strip tests and particle counting. *Clin Chem Lab Med* 2008; 46:703-13.
14. Agenzia Sanitaria e Sociale Regionale. Infezioni delle vie urinarie nell'adulto. Linea Guida Regionale. 190 – 2010. [www//asr.regione.emilia.romagna.it/wcm/asr/collana_dossier/doss190.htm](http://www.asr.regione.emilia.romagna.it/wcm/asr/collana_dossier/doss190.htm).
15. Bertolissi A, Lanzafame P, Screm M, Trevisan R, Rossi P, Bierti R, et al. Azienda Ospedaliera Santa Maria della Misericordia. Linee Guida: La Diagnosi Microbiologica delle infezioni Urinarie. <http://www.ospedaleudine.it>.
16. Gillespie T, Fewster J, Masterton R. The effect of specimen processing delay on borate urine preservation. *J Clin Pathol* 1999; 52:95-8.
17. Meers P, Chow C. Bacteriostatic and bactericidal actions of boric acid against bacteria and fungi commonly found in urine. *J Clin Pathol* 1990; 43:484-7.
18. Watson P, Duerden B. Laboratory assessment of physical and chemical methods of preserving urine specimens. *J Clin Pathol* 1977; 30:532-6.
19. Laurer B, Reller B, Mirrett S. Evaluation of preservative fluid for urine collected for culture. *J Clin Microbiol* 1979; 10:42-5.
20. Fromm P, Bieganiec B, Ehenrich Z, Barak M. Stability of common analytes in urine refrigerated for 24 h before automated analysis by test strips. *Clin Chem* 2000; 46:1384-6.

21. Eriksson I, Lindmann R, Thore M, Microbiological evaluation of a commercial transport system for urine samples. *Scand J Clin Lab Invest* 2002; 62:325-35.
22. Keuzer M, Brandts R. Flow cytometer and the urine Laboratory: Field evaluation of the Sysmex UF-100. *Sysmex Intern Journ* 1997; 7:117-22.
23. Ottiger C, Huber A. Quantitative urine particles analysis: Integrative approach for the optimal combination of automation with UF-100 and microscopic review with KOVA cell chambre. *Clin Chem* 2003; 49:617-23.
24. Lamchiagdhas P, Preechaborisutkull K, Lomsomboon P, Srisuchart P, Khan-u-Ra N, Preechaborisutkull B, et al. Urine sediment examination method and the iQ200 automated urine microscopy analyzer. *Clin Chim Acta* 2005; 358:167-74.
25. Wah D, Porntip K, Butch A. Analytic performance of the iQ200 automated urine microscopy analyzer and comparison with manual counts using Fuchs-Rosenthal cell chambers. *Am J Clin Pathol* 2005; 123:290-6.
26. Zaman Z, Fogazzi GB, Garigali G, Croci MD, Bayer G, Kránicz T. Urine sediment analysis: Analytical and diagnostic performance of sediMAX - a new automated microscopy image-based urine sediment analyser. *Clin Chim Acta* 2010; 41:147-54.
27. European Committee for Standardization. In vitro diagnostic system: transport packages for medical and biological specimens. Requirements tests: Brussels 1996.
28. Manoni F, Valverde S, Caleffi A, Alessio M, Silvestri M, de Rosa R. et al. Stability of common analytes and urine particles stored at room temperature before automated analysis. *Riv It Med Lab* 2008; 4:192-98.
29. Manoni F, Valverde S, Caleffi A, Alessio M, Gessoni G, Cappelletti P. Esame delle urine: è tempo di cambiare. *Riv It Med Lab* 2010; 6:52-8.
30. Tencer J, Thysell H, Andersson K, Grubb A. Stability of albumin, protein HC, immunoglobulin G, Kappa and Lambda chain immunoreactivity, arosomuroid and alpha 1 antitrypsin in urine stored at various condition. *Scand J Clin Lab Invest* 1994; 54:199-206.
31. Kierkegaard H, Feldt-Rasmussen U, Hoerder M, Andersen HJ, Jorgensen PJ. Falsely negative urinary leukocyte counts due to delayed examination examination. *Scand J Clin Lab Invest* 1980; 40:259-61.
32. Ampalahan M, Birch D, Becker G. Chimica preservation of urine sediment for phase contrast microscopy examination. *Nephron* 1994; 68:180-3.
33. Porter J, Brodie J. Boric acid preservation of urine samples. *Br Med J* 1969; 2: 353-5.
34. Grosso S, Bruschetta G, De Rosa R, Avolio M, Camporese A. Improving the efficiency of pre analytical and analytical work flow of urine cultures with flow cytometry. *New Microb* 2008; 31:501-5.
35. Roggeman S, Zaman Z. Safely reducing manual urine microscopy analyses by combining urine flow cytometer and strip results. *Am J Clin Pathol* 2001; 116(6):872-8.
36. Fuller C, Threatte G, Henry J. Basic examination of the urine. In: Henry JB, Davey FR, Herman CJ, et al., eds. *Clinical diagnosis and management by laboratory methods*. 20th ed. Philadelphia: WB Saunders, 2001:367-402.
37. Miller J. *Guide to specimen management in clinical microbiology*. American Society for Microbiology 1996.
38. Addis T. The number of formed elements in the urinary sediment of normal individuals. *J Clin Invest* 1926; 2:409-15.
39. Gadeholt H. Quantitative estimation of urinary sediment with special regard to source of errors. *Br Med J* 1964; 1:1547-9.
40. Fairley K, Birch D. Microscopic urinalysis in glomerulonephritis. *Kidney Int* 1993; 44:S9-12.
41. Shayanfar N, Tobler U, von Eckardstein A, Bestmann L. Automated urinalysis: first experiences and a comparison between the Iris iQ200 urine microscopy system, the

- Sysmex UF-100 flow cytometer and manual microscopic particle counting. *Clin Chem Lab Med* 2007;45:1251-6. Erratum in: *Clin Chem Lab Med* 2007; 45:1570.
42. Regeniter A, Haenni V, Risch L, Kochli HP, Colombo JP, Frei R, et al. Urine analysis performed by flow cytometry: reference range determination and comparison to morphological findings, dipstick chemistry and bacterial culture results--a multicenter study. *Clin Nephrol* 2001; 55(5):384-92.
 43. Kass E. Asymptomatic infections of the urinary tract. *Trans Assoc Am Phys* 1956; 69: 56-63.
 44. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract; with observations on the use of methionine as a urinary antiseptic. *AMA Arch Intern Med* 1957; 100:709-14.
 45. Lohr J, Donowitz L, Dudley S. Bacterial contamination rates for non-clean-catch and clean-catch midstream urine collections in boys. *J Pediatr* 1986; 109:659-60.
 46. Lohr JA, Donowitz LG, Dudley SM. Bacterial contamination rates in voided urine collections in girls. *J Pediatr* 1989; 114:91-3.
 47. Lifshitz E, Kramer L. Outpatient urine culture: does collection technique matter? *Arch Intern Med* 2000; 160:2537-40.
 48. Prandoni D, Boone M, Larson E, Blane C, Fitzpatrick H. Assessment of urine collection technique for microbial culture. *AmJ Infect Control* 1996; 24:219-21.
 49. Morris R, Watts M, Reeves D. Perineal cleansing before midstream urine, a necessary ritual? *Lancet* 1979; 2:158-9.
 50. Leisure M, Dudley S, Donowitz L. Does a clean-catch urine sample reduce bacterial contamination? *N Engl J Med* 1993; 328:289-90.
 51. Wheldon DB, Slack M. Multiplication of contaminant bacteria in urine and interpretation of delayed culture. *J Clin Pathol* 1977; 30:615-9.
 52. Rogers J, Saunders C. Urine collection in infants and children. *Nurs Times* 2008; 104(5):40-2.
 53. Koneman E, Allen S, Dowell V, Sommers H. *Testo Atlante di Microbiologia Diagnostica*, Antonio Delfino Roma 1987.
 54. Cheng Y, Wong S. Diagnosing symptomatic urinary tract infections in infants by catheter urine culture. *J Paediatr Child Health* 2005; 41:437-40.
 55. Gilljam B, Svensson M. In-out catheterization of young children with suspected urinary tract infection: a retrospective journal study. *Pediatr Nurs* 2008; 34(3):241-5.
 56. Metchock BG, Nolte FS, Wallace RJ. Mycobacterium. In: Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. *Manual of clinical microbiology*. 7th ed. Washington, DC: American Society for Microbiology Press. 1999:399-437.
 57. Peter J, Green C, Hoelscher M, Mwaba P, Zumla A, Dheda K. Urine for the diagnosis of tuberculosis: current approaches, clinical applicability, and new developments. *Curr Opin Pulm Med* 2010; 16:262-70.
 58. Anosike J, Okere A, Nwoke BE, Chukwu JU, Nwosu DC, Njoku-Tony RF, et al. Endemicity of vesical schistosomiasis in the Ebonyi Benue river valley, south eastern Nigeria. *Int J Hyg Environ Health* 2003; 206:205-10.

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